

Central Proteomics Facility

Tryptic Digest Protocol

The facility can perform digests for you (at a charge), or can provide the chemicals, equipment and facilities for you to perform them in the facility.

Solutions required:

Solution B: 0.04 g ammonium bicarbonate (ammonium hydrogen carbonate) in 10 mL Milli-Q grade water and 10 mL HPLC grade acetonitrile.

Acetonitrile: HPLC grade.

25 mM ammonium bicarbonate: 0.04 g in 20 mL Milli-Q grade water.

10 mM DTT: 0.031 g in 20 mL in 25mM ammonium bicarbonate solution.

55 mM Iodoacetamide: 0.2 g in 20 mL ammonium bicarbonate solution.

Promega Sequencing Grade Modified Trypsin (catalogue number: V5111. How to reconstitute the trypsin: to 1 vial add 100 uL Promega resuspension buffer. Store in 5 uL aliquots at -20 C until use. Prior to use add 95 uL 25 mM ammonium bicarbonate – the trypsin is now active. Any trypsin remaining after use should be disposed of.

Extraction buffer: 10 mL Milli-Q grade water + 10 mL acetonitrile + 20 microlitres formic acid.

Protocol

- 1) Cut-out gel bands on a clean surface with a clean scalpel. Dice into 1mm cubes and place in tube.
- 2) Wash bands with 100 uL Solution for 30 mins then discard supernatant

- 3) Repeat Wash with 100 uL Solution B 30 mins then discard supernatant (repeat until bands are no longer a strong blue colour).
- 4) Wash in 100 uL 100% acetonitrile for 10 mins and remove (gel pieces should dehydrate and go white).
- 5) Dry in SpeedVac for 10 mins or remove cap for 10 mins (all acetonitrile should have evaporated).
- 6) Add 100 uL 10 mM DTT for 30 minutes at 37 C and discard.
- 7) Wash gel pieces with 25mM ammonium bicarbonate solution and discard.
- 8) Wash gel pieces with 25mM ammonium bicarbonate solution and discard.
- 9) Wash in acetonitrile until gel pieces are white and discard.
- 10) Add 100 uL 55 mM Iodoacetamide for 60 mins (IN DARK) and discard
- 11) Wash bands with 100 uL Solution B for 10 mins and discard.
- 12) Wash bands with 100 uL Solution B for 10 mins and discard.
- 10) Wash in 100 uL 100% acetonitrile until gel pieces are white, then discard.
- 11) Dry in SpeedVac for 10 mins or remove cap for 10 mins.
- 12) Add 20 uL trypsin to each tube.
- 13) Digest overnight at 37 C.
- 14) Add 1 uL of formic acid to stop digest.
- 15) Remove supernatant containing peptides to clean new tube.
- 16) Add 50uL of extraction buffer to gel piece and incubate for 30 mins.
- 17) Remove supernatant containing peptides and pool with existing supernatant.
- 18) SpeedVac to dryness (available in Central Proteomics Facility) and bring to Central Proteomics Facility ASAP.

Central Proteomics Facility, Sir William Dunn Pathology School, South Parks Rd., Oxford,
OX1 3RE,

Tel: 01865 275 613, Email: Benjamin.thomas@path.ox.ac.uk